




NEWCLIP-TECHNICS

INNOVATION MEANS MOTION



ACTIV ANKLE
DISTAL AND DIAPHYSEAL FIBULA
MEDIAL MALLEOLUS

POLYAXIAL LOCKING FIXATION
DUALTEC SYSTEM® II

- ▶ Precontoured implants
- ▶ Polyaxiality of 20°
- ▶ 2 surgical approaches: lateral and posterolateral
- ▶ Medial malleolus fixation

ACTIV ANKLE

Indications: The Activ Ankle range is intended for the fixation of fractures, osteotomies and pseudarthroses of the distal and the diaphyseal fibula, the distal tibia and for the syndesmotic repair in adults.

Contra-indications:

- Serious vascular deterioration, bone devitalization.
- Pregnancy.
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone quality preventing the correct insertion of the implants into the bone.
- Muscular deficit, neurological deficiency or behavioral disorders, which could submit the implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

DISTAL FIBULA PLATES

LATERAL APPROACH



STANDARD PLATE

- Fixation of osteoporotic bones and complex fractures with or without syndesmosis injuries.



NARROW PLATE

- Fixation of simple fractures with or without syndesmosis injuries.

POSTEROLATERAL APPROACH



POSTEROLATERAL PLATE

- Fixation of short oblique fractures (Type B, as defined by the AO and Weber Classifications).

TECHNICAL FEATURES

A COMPREHENSIVE RANGE OF PLATES


LATERAL STANDARD

LATERAL NARROW


POSTEROLATERAL

DIAPHYSEAL PLATE

For diaphyseal simple fractures




Bendable sections




- Ⓢ : Ø3.5 mm cortical screws (CT3.5Lxx) or Ø4.0 mm lag screws (QT4.0Lxx) for the syndesmosis fixation.


- DTS2® polyaxial holes for locking and non-locking Ø2.8 mm screws.



- Holes for locking and non-locking Ø3.5 mm screws.



- Oblong holes for Ø3.5 mm cortical screws.



TECHNICAL FEATURES

A PRECONTOURED IMPLANT

→ OPTIMAL ANATOMICAL CONGRUENCE

The design of this implant is the result of a proprietary state-of-the-art mapping technology to establish the optimized congruence between the plate and the bone.




PLATE BENDING

The diaphyseal part of the fibula can vary from one patient to the other. In the case of long plates, it is possible to bend the diaphyseal part of plates at each appropriate area using bending irons (ANC542) for an optimal fit to the bone anatomy, following the instructions below:

- > Bending is only possible in the areas intended for this purpose,
- > A bendable area should be bent only once and in one direction,
- > Bending should not be performed excessively,
- > The holes must be protected so as to avoid damaging the fixation.





ANGULAR RANGE: +/- 10° POLYAXIAL LOCKING FIXATION

→ THE DTS2® TECHNOLOGY ALLOWS THE SCREW TO LOCK INTO THE PLATE WHILE PERMITTING AN ANGULATION OF THE SCREW.

Newclip Technics plates combine both polyaxial and locking technologies to create a fixed-angle construct particularly useful for poor bone quality and/or multifragmentary fractures.


The DTS2® polyaxial locking holes are located in the epiphyseal area, thus facilitating the insertion of the epiphyseal screws in diverging or converging directions and allowing for optimal strength of the assembly.



MONOAXIAL LOCKING FIXATION

→ FEATURES

- > The threaded sections under the screw head and inside the hole have **the same characteristics** (1):
 - Cylindrical internal thread profile,
 - Cylindrical external thread profile,
- > Screw head cap (2),
- > Plate and screw made from the same material: titanium alloy.




→ RESULTS

- > The screw is stopped in the hole by its cap, insuring the locking (3),
- > A perfect coaptation of both profiles during locking (4).

SURGICAL TECHNIQUE

PLACEMENT OF THE LATERAL PLATE

→ PRELIMINARY REDUCTION OF THE FRACTURE WITH A SCREW



Remark:

As an osteosynthesis screw used alone cannot bear weight and resist shear stresses, a plate should be used to allow early mobilization.


1. Reduce and temporarily maintain the fracture with bone reduction forceps (ANC504), making sure not to hinder the subsequent positioning of the screw. Drill with the Ø2.7 mm drill bit (ANC089C) using the guide gauge¹ (ANC191). The drilling should be perpendicular to the line of fracture.

2. When a lag effect is desired, over-drill the anterior cortex only by using the Ø3.5 mm drill bit (ANC542) according to the allowed compression principles. To simplify the procedure, it is also possible to use the reamer of the 2 in 1 instrument (ANC083C).

3. Insert the Ø3.5 mm cortical screw (CT3.5Lxx) through the line of fracture using the screwdriver part of the 2 in 1 instrument (ANC083C). In the case of osteoporotic bone, a compression washer (WASH-T4) can be added under the screw head so as to obtain optimized compression.

(1) The screw length can be directly read on the guide gauge. Always ensure that the guide gauge sits flush against the bone surface.

→ PLACEMENT OF THE PLATE




4. Check the positioning of the fast guide thanks to the 'DISTAL' and 'PROXIMAL' marks. Lock the fast guide onto the plate with the screwdriver (ANC082E).

5. The plate can be temporarily held in position with Ø1.3 mm pins (33.0213.120).

6. Drill (ANC089C) using the guide gauge (ANC191). The screw length can be directly read on the guide gauge.

7. Insert a Ø3.5 mm cortical screw (CT3.5Lxx) into the oblong hole with the screwdriver part of the 2 in 1 instrument (ANC083C). For optimal positioning, slide the plate using the oblong hole and tighten the cortical screw.



8. For the epiphyseal fixation, drill using the threaded guide gauge (ANC268C) for polyaxial fixation (a) or the non-threaded guide gauge (ANC046C) for monoaxial fixation (b) through the pre-angled fast guide. The screw length can be directly read on the guide gauge. For the monoaxial distal hole, drill (ANC088C) using the threaded guide gauge (ANC268C).

9. Insert a Ø2.8 mm locking screw (SDT2.8Lxx) through the fast guide using the screwdriver (ANC082E).

10. Repeat the whole procedure to insert the remaining distal Ø2.8 mm locking screws (SDT2.8Lxx) and remove the fast guide.

SURGICAL TECHNIQUE



11. For the diaphyseal fixation, drill using the guide gauge (ANC186) and insert the Ø3.5 mm locking screws (SOT3.5Lxx). For the Ø3.5 mm cortical screws (CT3.5Lxx) insertion, repeat this procedure using the guide gauge (ANC191)². The final tightening of the screws must be performed by hand.

NB: To ease the insertion of the Ø3.5 mm locking screws (SOT3.5Lxx), use the reamer part of the 2 in 1 instrument (ANC083C) to widen the first cortex previously drilled.

Remark:
The fixation steps remain unchanged for Narrow (RTSL-Nx) or Posterolateral plates (RTxQ1).

FINAL RESULT

(2) In the case of a bicortical fixation, the drilling depth can be checked on the length gauge (ANC124).

SYNDESMOSIS FIXATION




1. Drill (ANC256M) through the holes designed for syndesmosis screws using the guide gauge (ANC261M). The screw length can be directly read on the guide gauge.

2. Insert (ANC083C) a syndesmosis screw: Ø3.5 mm solid cortical screw (CT3.5Lxx) or Ø4.0 mm solid lag screw (QT4.0Lxx) into the appropriate oblong hole and/or standard hole designed for that purpose.

FINAL RESULT

MEDIAL MALLEOLUS FIXATION


→ **OPTION 1: WITH A CANNULATED SCREW**
(Ø4.0 mm cannulated compressive screw)



1. Insert the guiding Ø1.3 mm pin (33.0213.120). Then, introduce the Ø2.9mm cannulated drill bit (ANC414M) onto the guiding pin and drill. Read the screw length on the drill bit.

2. Insert the Ø4.0 mm compressive cannulated screw (H1.4QT4.0Lxx) using the cannulated screwdriver (ANC388) then remove the pin.

→ **OPTION 2: WITH SOLID SCREW**
(Ø3.5 mm solid cortical screw / Ø4.0 mm solid lag screw)



1. Drill (ANC089C) using the guide gauge (ANC191). The screw length can be directly read on the guide gauge. Always ensure that the guide gauge sits flush against the bone surface.


2. Insert the Ø4.0 mm lag screw (QT4.0Lxx) or the Ø3.5 mm cortical screw (CT3.5Lxx) using the screwdriver part of the 2 in 1 instrument (ANC083C).

Remarks:

1. In the case of osteoporotic bone, add a compression washer (WASH-T4) under the screw head so as to obtain optimal compression (see above).
2. Follow the whole procedure for adding a second screw.

IMPLANTS REFERENCES

→ DISTAL PLATES




LATERAL STANDARD PLATES

| Ref. | Description |
|--------|---|
| RTGLS1 | Lateral plate for distal fibula - Standard Left - Size 1 |
| RTDLS1 | Lateral plate for distal fibula - Standard Right - Size 1 |
| RTGLS2 | Lateral plate for distal fibula - Standard Left - Size 2 |
| RTDLS2 | Lateral plate for distal fibula - Standard Right - Size 2 |
| RTGLS3 | Lateral plate for distal fibula - Standard Left - Size 3 |
| RTDLS3 | Lateral plate for distal fibula - Standard Right - Size 3 |
| RTGLS4 | Lateral plate for distal fibula - Standard Left - Size 4 |
| RTDLS4 | Lateral plate for distal fibula - Standard Right - Size 4 |
| RTGLS5 | Lateral plate for distal fibula - Standard Left - Size 5 |
| RTDLS5 | Lateral plate for distal fibula - Standard Right - Size 5 |



LATERAL NARROW PLATES

| Ref. | Description |
|--------|---|
| RTSLN1 | Lateral plate for distal fibula - Narrow symmetrical - Size 1 |
| RTSLN2 | Lateral plate for distal fibula - Narrow symmetrical - Size 2 |



POSTEROLATERAL PLATES

| Ref. | Description |
|-------|---|
| RTGQ1 | Posterolateral plate for distal fibula - Left - Size 1 |
| RTDQ1 | Posterolateral plate for distal fibula - Right - Size 1 |
| RTGQ2 | Posterolateral plate for distal fibula - Left - Size 2 |
| RTDQ2 | Posterolateral plate for distal fibula - Right - Size 2 |
| RTGQ3 | Posterolateral plate for distal fibula - Left - Size 3 |
| RTDQ3 | Posterolateral plate for distal fibula - Right - Size 3 |

→ DIAPHYSEAL PLATE




DIAPHYSEAL PLATE

| Ref. | Description |
|------|---|
| FTS1 | Plate for diaphyseal fibula fracture - Size 1 |

IMPLANTS REFERENCES


→ Ø2.8 MM SCREWS



LOCKING SCREWS*

| Ref. | Description |
|-----------|--------------------------------|
| SDT2.8L10 | Ø2.8 mm locking screw - L10 mm |
| SDT2.8L12 | Ø2.8 mm locking screw - L12 mm |
| SDT2.8L14 | Ø2.8 mm locking screw - L14 mm |
| SDT2.8L16 | Ø2.8 mm locking screw - L16 mm |
| SDT2.8L18 | Ø2.8 mm locking screw - L18 mm |
| SDT2.8L20 | Ø2.8 mm locking screw - L20 mm |
| SDT2.8L22 | Ø2.8 mm locking screw - L22 mm |
| SDT2.8L24 | Ø2.8 mm locking screw - L24 mm |

* Green anodized.



NON LOCKING SCREWS*

| Ref. | Description |
|-----------|------------------------------------|
| QDT2.8L10 | Ø2.8 mm non locking screw - L10 mm |
| QDT2.8L12 | Ø2.8 mm non locking screw - L12 mm |
| QDT2.8L14 | Ø2.8 mm non locking screw - L14 mm |
| QDT2.8L16 | Ø2.8 mm non locking screw - L16 mm |
| QDT2.8L18 | Ø2.8 mm non locking screw - L18 mm |
| QDT2.8L20 | Ø2.8 mm non locking screw - L20 mm |
| QDT2.8L22 | Ø2.8 mm non locking screw - L22 mm |
| QDT2.8L24 | Ø2.8 mm non locking screw - L24 mm |

* Golden anodized.


Remark:

Please note that all implants are also available in sterile packaging.

An '-ST' code is added at the end of the reference.

Eg. : « SDT2.8L10-ST »


→ Ø3.5 MM SCREWS



LOCKING SCREWS*

| Ref. | Description |
|-----------|--------------------------------|
| SOT3.5L10 | Ø3.5 mm locking screw - L10 mm |
| SOT3.5L12 | Ø3.5 mm locking screw - L12 mm |
| SOT3.5L14 | Ø3.5 mm locking screw - L14 mm |
| SOT3.5L16 | Ø3.5 mm locking screw - L16 mm |
| SOT3.5L18 | Ø3.5 mm locking screw - L18 mm |
| SOT3.5L20 | Ø3.5 mm locking screw - L20 mm |
| SOT3.5L22 | Ø3.5 mm locking screw - L22 mm |
| SOT3.5L24 | Ø3.5 mm locking screw - L24 mm |


* Blue anodized.



NON-LOCKING SCREWS*

| Ref. | Description |
|-----------|------------------------------------|
| QOT3.5L10 | Ø3.5 mm non locking screw - L10 mm |
| QOT3.5L12 | Ø3.5 mm non locking screw - L12 mm |
| QOT3.5L14 | Ø3.5 mm non locking screw - L14 mm |
| QOT3.5L16 | Ø3.5 mm non locking screw - L16 mm |
| QOT3.5L18 | Ø3.5 mm non locking screw - L18 mm |
| QOT3.5L20 | Ø3.5 mm non locking screw - L20 mm |
| QOT3.5L22 | Ø3.5 mm non locking screw - L22 mm |
| QOT3.5L24 | Ø3.5 mm non locking screw - L24 mm |

* Purple anodized.




CORTICAL SCREWS*

| Ref. | Description |
|----------|---------------------------------|
| CT3.5L10 | Ø3.5 mm cortical screw - L10 mm |
| CT3.5L12 | Ø3.5 mm cortical screw - L12 mm |
| CT3.5L14 | Ø3.5 mm cortical screw - L14 mm |
| CT3.5L16 | Ø3.5 mm cortical screw - L16 mm |
| CT3.5L18 | Ø3.5 mm cortical screw - L18 mm |
| CT3.5L20 | Ø3.5 mm cortical screw - L20 mm |
| CT3.5L22 | Ø3.5 mm cortical screw - L22 mm |
| CT3.5L24 | Ø3.5 mm cortical screw - L24 mm |

* Not anodized or light blue anodized for sterile screws.


→ SYNDESMOSIS AND MEDIAL MALLEOLUS SCREWS



Ø3.5 mm CORTICAL SCREWS*

| Ref. | Description |
|----------|---------------------------------|
| CT3.5L40 | Ø3.5 mm cortical screw - L40 mm |
| CT3.5L45 | Ø3.5 mm cortical screw - L45 mm |
| CT3.5L50 | Ø3.5 mm cortical screw - L50 mm |
| CT3.5L55 | Ø3.5 mm cortical screw - L55 mm |
| CT3.5L60 | Ø3.5 mm cortical screw - L60 mm |
| CT3.5L65 | Ø3.5 mm cortical screw - L65 mm |
| CT3.5L70 | Ø3.5 mm cortical screw - L70 mm |


* Not anodized or light blue anodized for sterile screws.



Ø4.0 mm LAG SCREWS*

| Ref. | Description |
|----------|----------------------------|
| QT4.0L40 | Ø4.0 mm lag screw - L40 mm |
| QT4.0L45 | Ø4.0 mm lag screw - L45 mm |
| QT4.0L50 | Ø4.0 mm lag screw - L50 mm |
| QT4.0L55 | Ø4.0 mm lag screw - L55 mm |
| QT4.0L60 | Ø4.0 mm lag screw - L60 mm |
| QT4.0L65 | Ø4.0 mm lag screw - L65 mm |
| QT4.0L70 | Ø4.0 mm lag screw - L70 mm |

* Not anodized.



Ø4.0 mm CANNULATED SCREWS* (for medial malleolus only)**

| Ref. | Description |
|--------------|---|
| H1.4QT4.0L40 | Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - L40 mm |
| H1.4QT4.0L45 | Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - L45 mm |
| H1.4QT4.0L50 | Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - L50 mm |
| H1.4QT4.0L55 | Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - L55 mm |
| H1.4QT4.0L60 | Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - L60 mm |
| H1.4QT4.0L65 | Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - L65 mm |
| H1.4QT4.0L70 | Self-drilling compressive screw - Ø4.0 mm - cannulated Ø1.4 mm - L70 mm |

* Not anodized.

** Optional, as a replacement for QT4.0Lxx

→ COMPRESSION WASHER, OPTIONAL



WASH-T4: Washer

INSTRUMENTS REFERENCES

| INSTRUMENTS | | |
|-------------|--|-----|
| Ref. | Description | Qty |
| ANC046C | Ø2.0 mm non-threaded guide gauge for Ø2.8 mm screws | 1 |
| ANC082E | 2.0 mm quick coupling hexagonal prehensor screwdriver | 1 |
| ANC083C | 2-in-1: 2.5 mm hexagonal prehensor screwdriver - Ø3.5 mm countersink | 2 |
| ANC084 | Ø2.8 mm quick coupling countersink | 1 |
| ANC088C | Ø2.0 mm quick coupling drill bit - L125 mm | 2 |
| ANC089C | Ø2.7 mm quick coupling drill bit - L125 mm | 2 |
| ANC102 | Length gauge for Ø2.8 mm screws | 1 |
| ANC103 | 2.0 mm hexagonal non prehensor screwdriver | 1 |
| ANC107 | 2.5 mm quick coupling hexagonal non prehensor screwdriver | 1 |
| ANC124 | Length gauge for Ø3.5 mm cortical screws | 1 |
| ANC186 | Ø2.7 mm threaded guide gauge for Ø3.5 mm screws | 2 |
| ANC191 | Ø2.7 mm non-threaded bent guide gauge for Ø3.5 mm screws | 1 |
| ANC252 | Fast drilling guide for RTGLSx plates | 1 |
| ANC253 | Fast drilling guide for RTDLSx plates | 1 |
| ANC256M | Ø2.7 mm quick coupling drill bit - L180 mm | 1 |
| ANC261M | Ø2.7 mm non-threaded bent long guide gauge for Ø3.5 and Ø4.0 mm screws | 1 |
| ANC268C | Ø2.0 mm threaded guide gauge for Ø2.8 mm screws | 2 |
| ANC349 | 15 cm verbrugge forceps | 2 |
| ANC350 | Ø4.5 mm AO quick coupling handle - Size 1 | 2 |

| INSTRUMENTS | | |
|-------------|--|-----|
| Ref. | Description | Qty |
| ANC452 | Bending iron | 2 |
| ANC454 | Fast drilling guide for RTGQx plates | 1 |
| ANC455 | Fast drilling guide for RTDQx plates | 1 |
| ANC456 | Fast drilling guide for RTSLNx plates | 1 |
| ANC463 | Ø3.5 mm quick coupling countersink | 1 |
| ANC503 | 150 mm Reduction forceps | 1 |
| ANC504 | 150 mm pointed reduction forceps | 1 |
| ANC542 | Ø3.5 mm quick coupling drill bit - L125 mm | * |
| 33.0213.120 | Pin Ø1.3 L120 mm | 6 |
| A10407M | 12 cm pin for washers | * |
| 30920 | Prehensive plier | * |

*: Optional


| INSTRUMENTS FOR CANNULATED SCREWS (optional) | | |
|--|---|-----|
| Ref | Description | Qty |
| ANC388 | 2.5 mm quick coupling hexagonal non prehensor screwdriver - cannula Ø1.4 mm | 1 |
| ANC414M | Quick coupling drill bit Ø2.9 mm - cannula 1.4 mm - L125 mm | 1 |

REMOVAL SET

If you have to remove Activ Ankle implants, make sure to order the Newclip Technics removal set which includes the following instruments :


- ANC103 for Ø2.8 mm screws: 2.0 mm hexagonal non prehensor screwdriver
- ANC107 or ANC016 for Ø3.5 mm and Ø4.0 mm screws: 2.5 mm hexagonal screwdrivers
- ANC350: Ø4.5 mm AO quick coupling handle - Size 1

→ SET DESCRIPTION




INSTRUMENTS TRAY (ANC254/I2)

SCREWS RACK (ANC254/R)




BASE (ANC254/B)

ANC254/I3 for CT3.5Lxx and QT4.0Lxx



OR

ANC254/I3 for CT3.5Lxx and H1.4QT4.0Lxx



IMPLANTS TRAY (ANC254/I1)

The information presented in this brochure is intended to demonstrate a NEWCLIP TECHNICS product. Always refer to the package insert, product label and/or user instructions before using any NEWCLIP TECHNICS product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your NEWCLIP TECHNICS representative if you have questions about the availability of NEWCLIP TECHNICS products in your area.

NEWCLIP TECHNICS (HQ)

PA de la Lande Saint Martin
45 rue des Garotières
44115 Haute-Goulaine, France
+33 (0)2 28 21 23 25
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP TECHNICS IBERIA

Calle Frederic Mompou, 4b
Sant Just Desvern
08960 Barcelona, España
+34 938 299 526
contact@newclipiberia.com
www.newcliptechnics.com

NEWCLIP TECHNICS USA

Newclip USA
642 Larkfield Center
Santa Rosa CA 95403, USA
+1 707 230 5078
customerservice@newclipusa.com
www.newclipusa.com

NEWCLIP TECHNICS GERMANY

Newclip GmbH
Pröllstraße 11, D-86157 Augsburg,
KKK Bldg. 502, 3-18-1 Asakusabashi
Deutschland
+49 (0)821 650 749 40
info@newclipgmbh.com
www.newclipgmbh.de

NEWCLIP TECHNICS JAPAN

Newclip Technics Japan K.K.
Newclip Technics Japan K.K.
KKK Bldg. 502, 3-18-1 Asakusabashi
Taito-Ku, Tokyo, 111-0053, Japan
+81 (0)3 58 25 49 81
www.newcliptechnics.com

NEWCLIP TECHNICS AUSTRALIA

Newclip Australia
3B/11 Donkin Street
West End 4101, Australia
+61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com