


ALIANS MIDSHAFT

> KIT DESCRIPTION

INSTRUMENTS + IMPLANTS



ANC107 2.5 mm quick coupling hexagonal non prehensor screwdriver


ANC083C 2 in 1: 2.5 mm hexagonal prehensor screwdriver
- Ø3.5 mm countersink

ANC089C Ø2.7 mm quick coupling drill bit - L125 mm (x2)

ANC350 Ø4.5 mm AO quick coupling handle - Size 1 (x2)

ANC124 Length gauge for Ø3.5 mm cortical screws

Alians Midshaft Plates:
Locking plates for forearm
ETSIS1 Size 1
ETSIS2 Size 2
ETSIS3 Size 3
ETSIS4 Size 4



ANC251 18 cm verbrugge forceps (x2)

ANC191 Ø2.7 mm non threaded bent guide
gauge for Ø3.5 mm screws


ANC186 Ø2.7 mm threaded guide gauge
for Ø3.5 mm screws (x3)

Screw rack

SOT3.5Lxx Ø3.5 mm locking screws (L10 to 26 mm)

QOT3.5Lxx Ø3.5 mm non locking screws (L10 to 26 mm)


CT3.5Lxx Ø3.5 mm cortical screws (L10 to 26 mm)



> KIT DESCRIPTION

IMPLANTS


PLATES	
Ref.	Description
ETSIS1	Locking plate for forearm - Size 1 - L80 mm
ETSIS2	Locking plate for forearm - Size 2 - L100 mm
ETSIS3	Locking plate for forearm - Size 3 - L130 mm
ETSIS4	Locking plate for forearm - Size 4 - L160 mm



Ø3.5 MM LOCKING SCREWS*	
Ref.	Description

SOT3.5Lxx Locking screw - Ø3.5 mm - L10 to 26 mm
(2 mm incrementation)


* Blue anodized



Ø3.5 MM NON LOCKING SCREWS*	
Ref.	Description

QOT3.5Lxx Non locking screw - Ø3.5 mm - L10 to 26 mm
(2 mm incrementation)

* Fuchsia anodized



Ø3.5 MM CORTICAL SCREWS*	
Ref.	Description

CT3.5Lxx Standard cortical screw - Ø3.5 mm - L10 to 26 mm
(2 mm incrementation)

* Not anodized



REMARK

All implants are available in a sterile version.
An «-ST» is added at the end of the reference.
Ex: «CT3.5L10-ST»



NEWCLIP TECHNICS
PA de la Lande Saint Martin, 45 rue des Garottières
44115 Haute Goulaine, France
+33 (0)2 28 21 23 25
orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP USA
642 Larkfield Center
Santa Rosa CA 95403, USA
+1 707 230 5078
customerservice@newclipusa.com
www.newclipusa.com

NEWCLIP GmbH
Pröllstraße 11
D-86157 Augsburg, Deutschland
+49 (0)821 650 749 40
info@newclipgmbh.com
www.newclipgmbh.de

NEWCLIP Australia
3B/11 Donkin Street
West End 4101, Australia
+61 (0)2 81 886 110
solutions@newclipaustralia.com
www.newcliptechnics.com

NEWCLIP TECHNICS Japan K.K
KKK Bldg 502, 3-18-1 Asakusabashi
Taïto-Ku, Tokyo, 111-0053, Japan
+81 (0)3 58 25 49 81
Fax: +81 (0)3 58 25 49 86
www.newcliptechnics.fr

NEWCLIP TECHNICS IBERIA
Calle Frederic Mompou, 4b
Sant Just Desvern
08960 Barcelona, España
+34 938 299 526
contact@newclipiberia.com
www.newcliptechnics.com

The information presented in this brochure is intended to demonstrate a NEWCLIP TECHNICS product. Always refer to the package insert, product label and/or user instructions before using any NEWCLIP TECHNICS product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your NEWCLIP TECHNICS representative if you have questions about the availability of NEWCLIP TECHNICS products in your area.